NMKL - NordVal International c/o Institute of Marine Research P.O. box 1870 Nordnes, N-5817 Bergen, Norway www.nmkl.org



NordVal International Certificate

Issued for: Compact Dry ETB Method for the

Enumeration of Enterobacteriaceae

NordVal No.: 034

First approval date: 1 December 2008

Renewal date: 1 December 2020

Valid until: 1 December 2022

Manufactured by: Supplied by:

Nissui Pharmaceutical Co.Ltd, HyServe GmbH & Co. KG,

3-23-9 Ueno, Hechenrainerstr 24,

Taito-ku, Tokyo, 110-8736 82449 Uffing, Japan Germany

Compact Dry ETB

The principle of this method is growth of dedicated bacteria on specific chromogen media.

The performance of this method has been compared to the reference method ISO 21528-2:2004: "Microbiology of foods and animal feeding stuffs. Horizontal method for the detection and enumeration of *Enterobacteriaceae* - part 2: Colony Count Method".

The validation studies have been conducted by Campden BRI, UK, according to ISO 16140-2:2016 and NordVal International Protocol 1.

NordVal International concludes that Compact Dry ETB provides equivalent results to ISO 21528-2:2004 for a broad range of foods. The production of Compact Dry ETB is certified according to ISO 9001 and ISO 13485.

Date: 1/12 2020

Yours sincerely

Hilde Skår Norli

Chair of NordVal International

Eystein Oveland

NMKL Secretary-General



PRINCIPLE OF THE METHOD

Compact Dry ETB is a ready-to-use selective plate containing glucose for the enumeration of *Enterobacteriaceae*. Pre-treat the samples according to ISO 6687 or NMKL 91. An aliquot of 1 ml of an appropriate dilution is plated onto Compact Dry ETB plate. The plate is incubated at $37 \pm 1^{\circ}$ C and colonies (red/purple) were counted after $24 \pm 2h$.

FIELD OF APPLICATION

The method has been tested on enumeration of *Enterobacteriaceae* in a broad range of foods.

HISTORY

In 2007, the method was validated according to the ISO 16140:2003. Every two years until 2018 the method has been renewed without any additional studies.

In 2018 a renewal study was performed to comply with the requirements for relative trueness and accuracy profile in the new standard ISO 16140-2:2016. As the design of the Interlaboratory study (ILS) is the same for the 2003 and 2016 versions of ISO16140, the data from the ILS data of 2007 are re-evaluated using the new statistical approach outlined in ISO16140-2:2016.

COMPARISON STUDY

Relative trueness study

The trueness study is a comparative study between results obtained by the reference method and the results of the alternative method. Different categories, types and items were tested as shown in Table 1 below.

Table 1 – Categories and types tested

Category	Types			
	Dry milk product e.g. milk powder, powder for milk based desserts, dried infant formula			
Heat processed milk and dairy products	Dairy products e.g. ice-cream, yogurts, cream, hard cheese, soft cheese, raw milk cheese			
	Pasteurised milk products e.g. skimmed, semi-skimmed, full fat and flavoured milks	5		
Fresh produce and fruits	Cut ready to eat fruit e.g. fruit mixes, fruit juices	5		
	Cut ready to eat vegetables e.g. Bagged pre-cut salads and shredded carrot, cabbage, vegetable juices			
	Leafy greens/Sprouts e.g. soy, mung, alfalfa,	5		
Pow poultry and mosts	Fresh poultry cuts e.g. turkey breast, turkey fillet	5		
Raw poultry and meats (Combined category raw/ RTC meats and poultry)	Fresh mince e.g. lamb, beef, pork	5		
	Processed ready to cook e.g. frozen patties, marinated kebabs, seasoned chicken breasts	5		
Ready to eat foods	Ready to eat poultry e.g. turkey fillet, chicken sausage, pate	5		
(Combined category RTE/RTRH meats and poultry)	Cooked fish products e.g. prawns, terrine, pate, smoked fish	5		
	Cooked meat e.g. ham, salami, pate, corned beef	5		
	Ready to re-heat refrigerated food e.g. cooked chilled foods, rice and pasta, products			
Multi component foods or meal components	Ready to re-heat food frozen e.g. fries, pizza	5		
	Composite foods with substantial raw ingredients e.gpasta salads, sandwiches, deli-salads	5		




75 samples were analysed, whereof 60 samples were artificially contaminated and 15 samples were naturally contaminated.

The relative trueness is illustrated by the use of a Bland-Altman plot, i.e. the difference (bias) between paired samples analysed with the reference method and the alternative method respectively, plotted against the mean values obtained by the reference method. In the plot, Upper and Lower limits are included as the bias \pm 2 times the standard deviation of the bias.

The Bland-Altman Plot in Figure 1, illustrates the difference obtained in the enumeration of *Enterobacteriaceae* in foods by the alternative and the reference method, respectively.

Figure 1 Bland-Altman Plot of the food categories tested



For 'All Categories' there are five in 75 values (7%) which lie outside the CLs. This is a little more than the expectation of less than one in 20 (5%). Of the five points outside of the CLs, the data covered 3 different food categories, and 4 different inoculated strains. Although there was a general slight negative bias to the data, only one data point was outside the lower CL and 4 were outside the upper CL.

NordVal International considers the relative trueness for satisfactory.

ACCURACY PROFILE

The accuracy profile study is a comparative study between the results obtained by the reference method and the results of the alternative method. Each item used were artificially contaminated obtaining three target levels; low (10² cfu/g), medium (10⁴ cfu/g) and high (10⁶ cfu/g). Five test portions of each level of each item were analysed, resulting in 150 samples. The tested categories, types, items and inoculated strains are provided in the Table 2.



Table 2 - Categories, types and food items

Category	Types	Strain	Item
Dairy products	Pasteurised	E. coli CRA 1476 from dried milk	Pasteurised cream
	dairy products	Enterobacter agglomerans CRA 5613 from milk powder	Cream cheese
Fruits and vegetables	Fresh produce	E.hermanii CRA 7477 from sesame seeds	Ready to cook Vegetable preparation
		Citrobacter amalonaticus CRA 7458 from beansprouts	Vegetable juice
Raw poultry and meats (Combined category raw/ RTC meats and poultry)	Fresh meat	Salmonella Brandenberg CRA 1070 from beef	Pork mince
		Proteus mirabilis CRA 1588 from poultry	Raw bacon
Ready to eat foods (Combined category RTE/RTRH meats and poultry)	Cooked fish products e.g. prawns	E.coli CRA 2003 from fish	Fresh prawns
		Klebsiella oxytoca ATCC 15926	Fish pate
Multi component foods	Composite foods with raw ingredients	Hafnia alvei CRA 400 from sandwich	Sandwiches
		E. adecarboxylata CRA 5501 from skimmed milk powder	Cooked chilled rice

The statistical results and the accuracy profiles are provided in the Figures 2 to 6.



Figure 2 Dairy Products




Figure 3 Fruit and vegetable products






Figure 4 Meat and poultry




Figure 5 Ready to eat foods







Figure 6 Multi component foods



The observed profiles are within the 0.5 log AL or the recalculated AL limit calculated according to ISO16140-2:2015 section 6.1.3.3.

All the accuracy profiles fulfil the performance criteria after the permitted recalculation, and the alternative method is accepted as being equivalent to the reference method.

THE SELECTIVITY OF THE METHOD (INCLUSIVITY/EXCLUSIVITY)

The selectivity study was performed according to ISO 16140-2:2016.

<u>Inclusivity</u> is the ability of an alternative method to detect the target analyte from a wide range of strains. In the original study 32 strains were studied. One of the 32 strains failed to grow on Compact Dry ETB. In the renewal study from 2016, 18 of the 23 strains tested were detected by both methods. Those not detected by either method were *Erwinia amylovora* 8037 and *Erwinia herbicola* 7057. Three strains were detected by the reference method but not by the alternative method- these were: *Serratia liquefaciens* 10670, *Rahnella aqualatis* NCIMB 13365 and *Yersinia intermedia* 380.

Exclusivity is the lack of interference from a relevant range of non-target strains of the alternative method. In the original study, 21 of the 23 strains did not grow on either methods. The two strains that did grow in VRBGA included a strain of *Aeromonas hydrophila* (strain 4111) which appeared typical in this medium, a strain of *Vibrio parahaemolyticus* (strain 15737) which grew but was atypical in appearance on the ETB medium and which produced typical colonies in VRBGA, although growth was poor. One strain of *Pasteurella bettyae* yielded typical colonies by both methods whereas tests with other *Pasteurella* strains, including an additional *P. bettyae* strain showed inhibition of these bacteria by both media. *Pasteurella* spp belong to the family Pasteurellaceae and not the Enterobacteriaceae, and both members of these families are capable of fermenting glucose, and although their optimum growth temperature is 37°C most are fastidious in their growth requirements. However, unlike members of the Enterobacteriaceae *Pasteurella* spp. are oxidase-positive.



In the 2018 study, of the 10 exclusivity strains tested, one strain was detected by both the alternate method and by the reference method, *A. sobria* CRA 8390.

CONCLUSION OF THE COMPARISON STUDY

The results of the method comparison study showed that the Compact Dry ETB provide equivalent results to the reference method ISO 21528-2:2004. The lowest validated level is 2.0 cfu/g.

INTERALBORATORY STUDY


The interlaboratory study was conducted in November 2007. Ten laboratories analysed samples of pasteurised milk artificially contaminated with defined numbers of *Esherichia coli* and *Enterobacter aerogenes* according to ISO 21528-2:2004 and Compact Dry ETB respectively.

The obtained results (log cfu/g) is given in Table 3, and illustrated by an Accuracy Profile in Figure 7.

Table 3 The interlaboratory study results in log cfu/g

	Reference method		method		Upper	Lower		
Level	Median	S _R	Median	S _R	Bias	Level	Level	\pm AL
1	2.57	0.12	2.25	0.20	-0.32	-0.04	-0.60	0.50
2	3.62	0.13	3.49	0.19	-0.13	0.14	-0.40	0.50
3	4.58	0.069	4.48	0.12	-0.1	0.07	-0.27	0.50

Figure 7 Accuracy Profile of the interlaboratory study



The lowest level has a negative bias, and thus the lower level is below –AL.

CONCLUSION

According to the comparison and the interlaboratory study no substantial differences were found between the Compact Dry ETB method and the reference method ISO 21528-2:2004 for the enumeration of *Enterobacteriaceae*.