

Vortex Mixer



True vortex mixing with innovative touch



VM-96T (Persian Blue)



VM-96A




Specification

Model		VM-96A	VM-96T
Control	Speed range (rpm)	Max. 3000	Max. 3000
	Speed control	Scale	Scale
	Operating mode	Touch, Run-L, Run-H	Touch, Run
	Motion	Orbital	Orbital
	Orbit (dia., mm / inch)	4 / 0.16	4 / 0.16
Material	Platform	Silicone	Silicone
	Body	Polypropylene	Polypropylene
	Motor type	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)
Dimensions	Platform (W x D, mm / inch)	76 x 76 / 3 x 3	76 x 76 / 3 x 3
	Exterior (W x D x H, mm / inch)	154 x 210 x 84 / 6.1 x 8.3 x 3.3	148 x 159 x 77 / 5.8 x 6.3 x 3
	Net weight (kg / lbs)	3.0 / 6.6	2.6 / 5.7
Protection class (DIN EN 60529)		IP 42	-
Electrical requirements		AC 100V to 240V, 50/60Hz, 3.3A Max.	
Cat. No.	KR plug	AAH360215K	AAH361215BK
	US plug	AAH360215U	AAH361215BU

Accessories for Vortex Mixer



Platform (pop-off)



Platform (standard)
Included as standard on equipment purchase



Microplate Tray



Cat. No.	Description
00VMS0000063	Platform (Pop-off)
00VMS0000058	Platform (standard)
00VMS0000008JT	Microplate tray
AAA36520	Tube holder (up to $\Phi 50$ mm)
AAA36521	Tube holder ($\Phi 10$ mm x 16ea)
AAA36522	Tube holder ($\Phi 12$ mm x 9ea)
AAA36523	Tube holder ($\Phi 15$ mm x 9ea)
AAA36524	Tube holder ($\Phi 20$ mm x 5ea)
AAA36525	Tube holder ($\Phi 25$ mm x 4ea)
AAA36526	Tube holder ($\Phi 8$ mm x 18ea & $\Phi 10$ mm x 28ea)

