

LC Vue™ Series

Application Research Labs, Cold Chain & Storage Facilities, Production & Chemical Handling Areas

- **Smart Link* Support** Seamless integration between Jeio Tech equipment via the LC Connected App.
- **Zero Data Loss** High-performance F-RAM ensures data backup during network outages.
- **Enterprise-grade Security Wi-Fi** Supports WPA Enterprise for strict corporate network standards.
- **Clear OLED Display** Remains operational and visible even in freezers or extreme cold conditions.
- **Calibration Service** KOLAS (ISO/IEC 17025) accreditation available. (Optional)

LC Vue-Insight



Standard: 1x PT100 Sensor included (Expandable up to 2 channels)
Range: -200 ~ 200 °C
Power: 12V DC Adapter & Batteries
UPS Function: Auto-switch to battery
Accessory: Warning Light (Optional)



Specification

※ Features and specifications are subject to change without prior notice.

Model	LC Vue-Insight	
Measurement	Measurement Parameters ch1: PT100 Sensor (Standard) ch2: PT100 Sensor (Optional)	Measurement Range -200 ~ 200 °C (Accuracy: ±0.5 °C)
Software	LC Connected App support	
External DC Power	5 ~ 30 V (Battery operation capable)	
Dimension (W x D x H, mm)	165 x 50 x 115	
Accessories	Warning Light (ch4) optional	
Cat. No.	BEA50655001	